



## ASR100B120B-T3A Datasheet

- 50 mm up to 305 mm travel per axis
- Up to 12  $\mu\text{m}$  accuracy; 2  $\mu\text{m}$  repeatability; 85 mm/s speed
- Easily adjustable non-contact hall sensor limit switches
- Supported by iManager and MetaMorph • Microscopy Automation & Image Analysis Software
- Designed for use with an X-MCC2 Series stepper motor controller or any 2-phase stepper motor controller
- With AutoDetect, the X-MCC controller configures its settings automatically for the connected peripheral
- Custom versions available

## ASR Series Overview

Zaber's ASR series microscope stages combined with our X-MCC2 series multi-axis controllers are designed as replacements for manual stages on upright and inverted microscopes or for stand-alone operation as scanning stages. Set up is easy with AutoDetect. Once connected, the X-MCC controller will automatically detect and configure the ASR.

Mounting adaptors are available for breadboards and most common microscopes. Custom adaptors and plates are available upon request.

*For more information visit: <https://www.zaber.com/products/scanning-microscope-stages/ASR>*

## ASR Series Part Numbering

ASR100B120B-T3A Drawings

## ASR100B120B-T3A Specifications

Specification	Value	Alternate Unit
Microstep Size (Default Resolution)	0.15625 $\frac{1}{2}$ $\frac{1}{2}$ m	
Built-in Controller	No	
Recommended Controller	X-MCC2 (48 V) Recommended	
AutoDetect	Yes	
Lower Travel Range	100 $\frac{1}{2}$ mm	3.937"
Upper Travel Range	120 $\frac{1}{2}$ mm	4.724"
Accuracy (unidirectional)	40 $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$ m	0.001575"
Repeatability	<1 $\frac{1}{2}$ 2 $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$ m	<1 $\frac{1}{2}$ 0.000079"
Flatness	20 $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$ m	0.0008"
Backlash	<1 $\frac{1}{2}$ 4 $\frac{1}{2}$ $\frac{1}{2}$ $\frac{1}{2}$ m	<1 $\frac{1}{2}$ 0.000157"
Maximum Speed	85 $\frac{1}{2}$ mm/s	3.346"/s
Minimum Speed	0.000095 $\frac{1}{2}$ mm/s	0.000004"/s
Speed Resolution	0.000095 $\frac{1}{2}$ mm/s	0.000004"/s
Peak Thrust	95 $\frac{1}{2}$ N	21.3 $\frac{1}{2}$ lb
Maximum Continuous Thrust	95 $\frac{1}{2}$ N	21.3 $\frac{1}{2}$ lb
Maximum Centered Load	100 $\frac{1}{2}$ N	22.4 $\frac{1}{2}$ lb
Guide Type	Crossed roller bearing	
Pitch	0.02 $\frac{1}{2}$	0.349 $\frac{1}{2}$ mrads
Roll	0.02 $\frac{1}{2}$	0.349 $\frac{1}{2}$ mrads
Yaw	0.005 $\frac{1}{2}$	0.087 $\frac{1}{2}$ mrads
Linear Motion Per Motor Rev	2.00 $\frac{1}{2}$ mm	0.079"
Motor Steps Per Rev	200	
Motor Type	Stepper (2 phase)	
Motor Rated Current	1500 $\frac{1}{2}$ mA/phase	
Motor Winding Resistance	2.05 $\frac{1}{2}$ ohms/phase	
Inductance	1 $\frac{1}{2}$ mH/phase	
Motor Connection	D-sub 15	
Default Resolution	1/64 $\frac{1}{2}$ of a step	
Mechanical Drive System	Precision lead screw	
Limit or Home Sensing	Magnetic, adjustable home and away sensors	

Specification	Value	Alternate Unit
Axes of Motion	2	
Mounting Interface	Separate mounting adaptors available	
Operating Temperature Range	0 to 50½½C	
RoHS Compliant	Yes	
CE Compliant	Yes	
Vacuum Compatible	No	
Weight	3½kg	6.614½lb

## ASR Series Charts





## Contact

Email: [contact@zaber.com](mailto:contact@zaber.com)

Phone (toll free Canada/USA): 1-888-276-8033

Phone (direct): 1-604-569-3780

Fax: 1-604-648-8033

Zaber Technologies Inc.

#2 - 605 West Kent Ave. N.

Vancouver, British Columbia

Canada, V6P 6T7

<https://www.zaber.com>